



# Recombinant *Xenopus laevis* SUMO-activating enzyme subunit 1 (sae1)

<b>Product Code</b>	CSB-EP809048XBE-B
<b>Abbreviation</b>	sae1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8JGT5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVEKEEAVIS EEEAAQYDRQ IRLWGLEAQK RLRTSRVLLV GMRGLGAEVA KNLILAGVKA LTLLDHEQVS SEDSRAQFLI PSGSLGQNRA EASLNRARNL NPMVSVEADT ENINQKSDDF FTQFDVVCLT SCPSDLLVRV NHICHKHNK FFTGDVYGYH GSMFADLGEH EFVEEKAKVT KAKPLVEDGP EAKKAKIDPT ETILVKKKVQ FCPLKDALEI DWRSEKAKSA LKKTPTDYFL LQVLMKFRTD KGRDPQPSSY QEDSELLLI CSDVLDLGLV SPDLLPKDFA SYCFSEMAPV CAVVGVLGQ EIVKALSLRD APHNNFFFFD GKTSNGIVDC LGSK
<b>Source</b>	<i>E. coli</i>
<b>Target Names</b>	sae1
<b>Protein Names</b>	Recommended name: SUMO-activating enzyme subunit 1 Alternative name(s): SUMO-activating enzyme E1 N subunit Ubiquitin-like 1-activating enzyme E1A
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.