



Recombinant Mouse Thiamine-triphosphatase (Thtpa)

Product Code	CSB-BP809067MO
Abbreviation	Thtpa
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8JZL3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AQGLIEVER KFAPGPDTEE RLQELGATLE HRVTFRDITYY DTSELSLMLS DHWLRQREGS GWELKCPGVT GVSOPHNEYV EVTSEAAIVA QLFELLGSGE QKPAGVAAVL GSLKLQEVAS FITTRSSWKL ALSGAHQEP QLTIDLDSAD FGYAVGEVEA MVHEKAEVPA ALEKIITVSS MLGVPAQEEA PAKLMVYLQR FRPLDYQRLLEAASSGEATG DSAS
Source	Baculovirus
Target Names	Thtpa
Protein Names	Recommended name: Thiamine-triphosphatase Short name= ThTPase EC= 3.6.1.28
Expression Region	2-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.