



# Recombinant Mouse Endonuclease 8-like 1 (Neil1)

<b>Product Code</b>	CSB-YP809118MO
<b>Abbreviation</b>	Neil1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K4Q6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PEGPELHLA SHFVNETCKG LVFGGCVEKS SVSRNPEVVPF ESSAYHISAL ARGKELRLTL SPLPGSQPPQ KPLSLVFRFG MSGSFQLVPA EALPRHAHLR FYTAPPAPRL ALCFVDIRRF GHWDPGGEWQ PGRGPCVLL EYERFRENVL R NLSDKAFDRP ICEALLDQRF FNGIGNYLRA EILYRLKIPP FEKARTVLEA LQQCRPSPEL TLSQKIKAKL QNPDLLELCH LVPKEVVQLG GKGYPGERGE EDFAAFRAWL RCYGVPGMSS LRDRHGRTIW FQGDGPPLAP KGGRSQKKKS QETQLGAEDR KEDLPLSSKS VSRMRRARKH PPKRIAQQSE GAGLQQNQET PTAPEK GKRR GQRASTGHRR RPKTIPDTRP REAGESSAS
<b>Source</b>	Yeast
<b>Target Names</b>	Neil1
<b>Protein Names</b>	Recommended name: Endonuclease 8-like 1 EC= 3.2.2.- EC= 4.2.99.18 Alternative name(s): DNA glycosylase/AP lyase Neil1 DNA-(apurinic or apyrimidinic site) lyase Neil1 Endonuclease VIII-like 1 Nei homolog 1 Short name=
<b>Expression Region</b>	2-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.