



Recombinant Isocitrate dehydrogenase kinase/phosphatase (aceK)

Product Code	CSB-EP809931EGX
Abbreviation	aceK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FB61
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MPRGLELLIA QTILQGFDAQ YGRFLEVTSQ AQQRFEQADW HAVQQAMKNR IHLYDHHVGL VVEQLRCITN GQSTDAAFL RVKEHYTRLL PDYPRFEIAE SFFNSVYCRL FDHRS LTPER LFIFSSQPER RFRTIPRPLA KDFHPDHGWE SLLMRVISDL PLRLRWQNKS RDIHYIVRHL TETLGTDNLA ESHLQVANEL FYRNKAAWL V GKLITPSGTL PFLPIHQTD DGELFIDTCL TTAEASIVF GFARSYFMVY APLPAALVEW LREILPGKTT AELYMAIGCQ KHAKTESYRE YL VYLQGCNE QFIEAPGIRG MVMLVFTLPG FDRVFKVIKD RFAPQKEMSA AHVRACYQLV KEHDRVGRMA DTQEFENFVL EKRHISPALM ELLLQEA AEK ITDLGEQIVI RHLYIERRMV PLNIWLEQVE GQQLRDAIEE YGN AIRQLAA ANIFPGDMLF KNFGVTRHGR VVFYDYDEIC YMTEVNFRDI PLPRYPEDEL ASEPWYSVSP GDVFPEEFRH WLCADPRIGP LFEEMHADLF RADYWRALQN RIREGHVEDV YAYRRRQRFS VRFV
Source	E.coli
Target Names	aceK
Protein Names	Recommended name: Isocitrate dehydrogenase kinase/phosphatase Short name= IDH kinase/phosphatase Short name= IDHK/P EC= 2.7.11.5 EC= 3.1.3.-
Expression Region	1-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.