



Recombinant Xaa-Pro dipeptidase (pepQ)

Product Code	CSB-MP809936EGX
Abbreviation	pepQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBI1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MESLASLYKN HIATLQERTR DALARFKLDA LLIHSGELFN VFLDDHPYPF KVNPQFKAWV PVTQVPNCWL LVDGVNPKPL WFYLPVDYWH NVEPLNSFW TEDVEVIALP KADGIGSLLP AARGNIGYIG PVPERALQLG IEASNINPKG VIDYLHYRS FKTEYELACM REAQKMAVNG HRAAEEAFRS GMSEFDINIA YLTATGHRDT DVPYSNIVAL NEHA AVLHYT KLDHQAPEEM RSFLLDAGAE YNGYAADLTR TWSAKSDNDY AQLVKDVNDE QLALIATMKA GVS YVDYHIQ FHQRIAKLLR KHQITDMSE EAMVENDLTG PFMPHGIGHP LGLQVHDVAG FMQDDSGTHL AAPAKYPYLR CTRLQPGMV LTIEPGIYFI ESLLAPWREG QFSKHFNWQK IEALKPFGGI RIEDNVVIHE NNVENMTRDL KLA
Source	Mammalian cell
Target Names	pepQ
Protein Names	Recommended name: Xaa-Pro dipeptidase Short name= X-Pro dipeptidase EC= 3.4.13.9 Alternative name(s): Imidodipeptidase Proline dipeptidase Short name= Prolidase
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.