



Recombinant 4-alpha-L-fucosyltransferase (wecF)

Product Code	CSB-YP809937EGX
Abbreviation	wecF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBP9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MTVLIHVLGS DIPHHNRTVL RFFNDALAAT SEHAREFMVA GKDDGLSDSC LALSQQFFPG KKSLAEAVIA KAKANRQQRFFHGGQFNPKL WLALLSGGIK PSQFFWHIWG ADLYELSSGL RYKLFYPLRR LAQKRVGCVF ATRGDLSFFA KTHPKVRGEL LYFPTRMDPS LNTMANDRQR EGKMTILVGN SGDRSNEHIA ALRAVHQQFG DTVKVVVPMG YPPNEAYIE EVRQAGLELF SEENLQVLSE KLEFDAYLTL LRQCDLGYFI FARQQGIGTL CLLIQAGIPC VLNRENPFWQ DMTEQHLPVL FTTDDLNEDI VREAQRQLAS VDKNTIAFFS PNYLQGQWQRA LAIAGEVA
Source	Yeast
Target Names	wecF
Protein Names	Recommended name: 4-alpha-L-fucosyltransferase EC= 2.4.1.- Alternative name(s): TDP-Fuc4NAc:lipid II Fuc4NAc transferase Short name= Fuc4NAc transferase
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.