



# Recombinant Exodeoxyribonuclease 7 large subunit (xseA)

<b>Product Code</b>	CSB-EP809957EGX
<b>Abbreviation</b>	xseA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FF64
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLPSQSPAIF TVSRLNQTVR LLEHEMGQV WISGEISNFT QPASGHWYFT LKDDTAQVRC AMFRNSNRRV TFRPQHGGQV LVRANITLYE PRGDYQIIVE SMQPAGEGLL QQKYEQLKAK LQAEGLFELQ YKKSPLSPAHCVGVITSKTG AALHDILHVL KRRDPSLPVI IYSTAVQGDD APGQIVRAIE LANKRNECDV LIVGRGGGSL EDLWSFNDER VARAIFASRI PIVSAVGHET DVTIADFVAD LRAPTPSAAA EVVSRNQQEL LRQVQSTRQR LEMAMDYYLA NRTRRFTQIH HRLQQQHPQL RLARQQTMLE RLQKRMSFAL ENQLKRAGQQ QQRLTRQLVQ QNPQSRIHRA QTRIQQLEYR LAETLRAQLS ATRERFGNAV THLEAVSPLS TLARGYSVTS AADGAVLKQV KQVKVGETLT TRLGDGVVIS EVSAVTKTRK SRKKT SNP
<b>Source</b>	E.coli
<b>Target Names</b>	xseA
<b>Protein Names</b>	Recommended name: Exodeoxyribonuclease 7 large subunit EC= 3.1.11.6 Alternative name(s): Exodeoxyribonuclease VII large subunit Short name= Exonuclease VII large subunit
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.