



Recombinant Alkanesulfonate monooxygenase (ssuD)

Product Code	CSB-EP809986EGX
Abbreviation	ssuD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FJ93
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	SLNMFWFLP THGDGHYLG T EEGSRPIDHG YLQQIAQAAD RLG YTGVLIP TGRSCEDAWL VAASMIPVTQ RLKFLVALRP SVTSPTVAAR QAATLDRLSN GRALFNLVTG SDPQELAGDG VFLDHSE RYE ASAEFTQVWR RLLLG ETVDF NGKHIHV RGA KLLFPPIQQP YPP LYFGGSS DVAQELAAEQ VDL YLTWGE P PELVKEKIEH VRAKAAAHGR KIRFGVRLHV IVRETND EAW QAAERLISRL DDETIAKAQA AFARTDSVGQ QRMAALHNGK RDNLEISP NL WAGVGLVRGG AGTALVGDGP TVAARINEYA ALGIDSFVLS GYPHLEEAYR VGELLFP HLD VAIPEIPQPQ PLNPQGEAVA NDFIPRNVAQ S
Source	E.coli
Target Names	ssuD
Protein Names	Recommended name: Alkanesulfonate monooxygenase EC= 1.14.14.5 Alternative name(s): FMNH ₂ -dependent aliphatic sulfonate monooxygenase
Expression Region	2-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.