



# Recombinant UvrABC system protein B (uvrB)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP809990EGX  |
| <b>Abbreviation</b>      | uvrB   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q8FJP8   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | SKPFLNSA FKPSGDQPEA IRRLEEGLED GLAHQTL LGV TGSGKTFTIA<br>NVIADLQRPT MVLAPNKTLA AQLYGEMKEF FPENAVEYFV SYYDYYQPEA<br>YVSSDTFIE KDASVNEHIE QMRLSATKAM LERRDVVVA SVSAIYGLGD<br>PDLYLKMMMLH LTVGMIIDQR AILRRLAELQ YARNDQAFQR GTFRVRGEVI<br>DIFPAESDDI ALRVXLFDDEE VERLSLFDPL TGQIVSTIPR FTIYPKTHYV<br>TPRERIVQAM EEIKEELAAR RKVLENNKL MEEQRLTQRT QFDLEMMNEL<br>GYCSGIENYS RFLSGRGPGE PPPTLFDYLP ADGLLVVDES HVTIPQIGGM<br>YRGDRARKET LVEYGFRLPS ALDNRPLKFE EFEALAPQTI YVSATPGNYE<br>LEKSGGDVVD QVVRPTGLLD PIIEV RPVAT QVDDLLSEIR QRAAINERVL<br>VTTLTKRMAE DLTEYLEEHG ERVRYLHSDI DTVERMEIIR DLRLGEFDVL<br>VGINLLREGL DMPEVSLVAI LDADKEGFLR SERSLIQTIG RAARNVNGKA<br>ILYGDKITPS MAKAIGETER RREKQQKYNE EHGITPQGLN KKVVDILALG<br>QNIAKTKAKG RGKSRPIVEP DNVPM DMSPK ALQQKIHELE GLMMQHAQNL<br>EFEEAAQIRD QLHQLRELF I AAS |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | uvrB   |
| <b>Protein Names</b>     | Recommended name: UvrABC system protein B Short name= Protein UvrB<br>Alternative name(s): Excinuclease ABC subunit B  |
| <b>Expression Region</b> | 2-673  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,  |



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.