



Recombinant Ribosomal RNA small subunit methyltransferase A (rsmA)

Product Code	CSB-EP810004EGX
Abbreviation	rsmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FL96
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNNRVHQGHL ARKRFGQNFL NDQFVIDSIV SAINPQKQQA MVEIGPGLAA LTEPVGERLD QLTVIELDRD LAARLQTHPF LGPKLTIYQQ DAMTFNFGEL AAKMGQPLRV FGNLPYNIST PLMFHLFSYT DAIADMHFML QKEVVNRLVA GPNSKAYGRL SVMAYYCNV IPVLEVPPSA FTPPPKVDSA VVRLVPHATM PHPVKDVRVL SRITTEAFNQ RRKTIRNSLG NLFVSEVLTG MGIDPAMRAE NISVAQYCQM ANYLAENAPL QES
Source	E.coli
Target Names	rsmA
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase A EC=2.1.1.182 Alternative name(s): 16S rRNA (adenine(1518)-N(6)/adenine(1519)-N(6))-dimethyltransferase 16S rRNA dimethyladenosine transferase 16S rRNA dimethylase
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.