



Recombinant Transaldolase B (talB)

Product Code	CSB-YP810006EGX
Abbreviation	talB
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FLD1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	<p>TDKLTSLRQ YTTVVADTGD IAAMKLYQPQ DATTNPSLIL NAAQIPEYRK LIDDAVAWAK QQSNDRAQQI VDATDKLAVN IGLEILKLVP GRISTEVDAR LSYDTEASIA KAKRLIKLYN DAGISNDRIL IKLASTWQGI RAAEQLEKEG INCNLTLIFS FAQARACAEA GVFLISPFVG RILDWYKANT DKKEYAPAED PGVSVSEIY QYYKEHGYET VVMGASFRNI GEILELAGCD RLTIPALLK ELAESEGAIE RKLSYTGEVK ARPARITese FLWQHNQDPM AVDKLAEGIR KFAVDQEKLE KMIGDLL</p>
Source	Yeast
Target Names	talB
Protein Names	Recommended name: Transaldolase B EC= 2.2.1.2
Expression Region	2-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>