



Recombinant Lipoyl synthase (lipA)

Product Code	CSB-YP810012CQG
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FNP4
Product Type	Recombinant Protein
Immunogen Species	Corynebacterium efficiens (strain DSM 44549 / YS-314 / AJ 12310 / JCM 11189 / NBRC 100395)
Purity	≥85% (SDS-PAGE)
Sequence	MTIAPEGRRM LRVEARNSQT PIETKPRWIR NQVKNNGPEYK DMKERVAGAS LHTVCQEAGC PNIHECWESR EATFLIGGAN CSRRCDFCMI NSARPEPLDR GEPLRVAESV REMRLNYSTI TGVTRDDLDPD EGAWLYSEVV RKIHELNPNT GVENLVPDFS GKRDLLQEVF ESRPEVFAHN LETVPRIFKR IRPAFRYDRS LDVIRQARDF GLVTKSNLIL GMGETREEIS EALRDLHSAG TDIITITQYL RPGPLFHPIE RWVKPEEFLE LSDEAYEIGF AAVMSGPLVR SSYRAGKLYA QALKFRGEEL PANLSHLAET TEGPTTQEAS SLLERYGASE DAPVIPRG
Source	Yeast
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.