



Recombinant Arabidopsis thaliana Probable E3 ubiquitin-protein ligase BAH1-like (RF178)

Product Code	CSB-MP810152DOA
Abbreviation	RF178
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GW10
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKFGETFTEY LHGEEWFLE KCRFVEYKKL KKVLLKCKTC NSTKSDDGQI IPSATSSSL S DSCECKACPW CDQMFFEELM KEATDIAGFF RSRVRHLLHL HVATGMQRYM IRLRRCFTDE KQALVQEGQI LIQYITMNAI AIRKILKKYD KVHSSSENGKN FKLKMRAERI ELLHSPWLIE LGAFYLN SGL DNVGNFKNSF GRVACENLNE DQPVLKMLP NSIELEYDLT CAICLET VFN PYALKCGHIF CNSCACSAAS VLIFQGIKAA PRHSKCPICR EAGVYAEAVH MIELHLL LKT RSKEYWKERM MNERSEMVKQ SKMFWNEQTK HMIGY
Source	Mammalian cell
Target Names	RF178
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase BAH1-like EC=6.3.2.- Alternative name(s): RING finger protein 178
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.