



Recombinant Arabidopsis thaliana GATA transcription factor 24 (GATA24)

Product Code	CSB-EP810162DOA
Abbreviation	GATA24
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GXL7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDDLHGRNGR MHIGVAQNPM HVQYEDHGLH HIDNENSMMMD DHADGGMDEG VETDIPSHPG NSADNRGEVV DRGIENGDQL TLSFQQQVYV FDRVSPEKVQ AVLLLLGGRE VPHTLPTTLG SPHQNNRVLG LSGTPQRLSV PQRLASLLRF REKRKGRNFD KTIRYTVRKE VALRMQRKKG QFTSAKSSND DSGSTGSDWG SNQSWAVEGT ETQKPEVLCR HCGTSEKSTP MMRRGPDGPR TLCNACGLMW ANKGTLRDLS KVPPPQTPQH LSLNKNEDAN LEADQMMEVT GDISNTQ
Source	E.coli
Target Names	GATA24
Protein Names	Recommended name: GATA transcription factor 24 Alternative name(s): Protein TIFY 2B ZIM-like 1 protein
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.