



# Recombinant Arabidopsis thaliana E3 SUMO-protein ligase MMS21 (MMS21)

<b>Product Code</b>	CSB-YP810169DOA
<b>Abbreviation</b>	MMS21
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8GYH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASASSSDGV AGRIQNASLV LVSDNSSTLA DIRKAVAMMK NIAVQLEKEN QTDKVKDLEN SVAELLDLHS DCNHRSTAIQ SVANRYQPVE QLTDFKLLD DEFTKLKATP SSVQNDHLM RQFREAVWNV HHAGEPMPGD DDEDIVMTST QCPLLNMTCP LSGKPVTELA DPVRSMDCRH VYEKSVILHY IVNNPNANCP VAGCRGKLQN SKVICDAMLK FEIEEMRSLN KQSNRAEVIE DFTEDVDED
<b>Source</b>	Yeast
<b>Target Names</b>	MMS21
<b>Protein Names</b>	Recommended name: E3 SUMO-protein ligase MMS21 EC= 6.3.2.- Alternative name(s): MMS21 homolog Short name= AtMMS21 Protein HIGH PLOIDY 2
<b>Expression Region</b>	1-249
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.