



Recombinant *Oryza sativa* subsp. japonica IAA-amino acid hydrolase ILR1-like 9 (ILL9)

Product Code	CSB-MP810186OFG
Abbreviation	ILL9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H3C7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	SAW AAAAGDDDLL AAAREPGMAE WLRGVRRRIH RHPELAFEEV RTSELVRAEL DAIGVPYQWP VARTGVVATI AGGGGGDGPV VALRADMDAL PVQELVDWEH KSQENGMHA CGHDAHTAML LGAAKLLQKR KNELKGTVKL VFQPAEEGSA GAYYVLQEGV LDDVSAMFGM HVDPALPVGVAARPGPFAA TSGRFLATIT GKGGHAAFPH DAIDPVVAAS NAILSLQQIV AREIDPLQGA VVSITFVKGG EAYNVIPQSV EFGGTMRSM DEGLAYLMKR IKEIVEGQAA VNRCGGGVDF MEESMRPYPA VVNDGMYAH ARASAERLLG AGGVRVAPQL MGAEDFGFYA ARMPSAFFTI GVGNATTSSA RAAHTTHSPH FVIDEALPV GAAVHAAVAI DYLSKHASSM
Source	Mammalian cell
Target Names	ILL9
Protein Names	Recommended name: IAA-amino acid hydrolase ILR1-like 9 EC= 3.5.1.-
Expression Region	28-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.