



Recombinant Human Keratin-associated protein 11-1 (KRTAP11-1)

Product Code	CSB-MP810266HU
Abbreviation	KRTAP11-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUC1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSFNCSTRNC SSRPIGGRCI VPVAQVTTTS TTDADCLGGI CLPSSFQTGS WLLDHCQETC CEPTACQPTC YRRTSCVSNP CQVTCRQTT CISNPCSTTY SRPLTFVSSG CQPLGGISSV CQPVGGISTV CQPVGGVSTV CQPACGVSRT YQQSCVSSCR RTC
Source	Mammalian cell
Target Names	KRTAP11-1
Protein Names	Recommended name: Keratin-associated protein 11-1 Alternative name(s): High sulfur keratin-associated protein 11.1
Expression Region	1-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.