



Recombinant Human Spermatogenesis-associated protein 2-like protein (SPATA2L)

Product Code	CSB-EP810268HU-B
Abbreviation	SPATA2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUW3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGSSSLSEEDY RQCLERELRR GRAGVCGDPS LRAVLWQILV EDFDLHGALQ DDALALLTDG LWGRADLAPA LRGLARAFEL LELAAVHLYL LPWRKEFTTI KTFSGGYVHV LKGVLSDDL LKSFQKMGYV RRDSHRLMVT ALPPACQLVQ VALGCFALRL ECEILGEVLA QLGTSVLP AE ELLQARRASG DVASCVAWLQ QRLAQDEEPP PLPPRGSPAA YRAPLDLYRD LQEDEGSEDA SLYGEPSPGP DSPPAELAYR PPLWEQSAKL WGTGGRAWEP PAEELPQASS PPYGALEEGL EPEPSAFSFL SLRRELSRPG DLATPESSAA ASPRRIRAEG VPASAYRSVS EPPGYQAHSC LSPGALPTLC CDTQRQLHAA HCAALPACRP GHSLRVLLGD AQRRLWLQRA QMDTLLYNSP GARP
Source	E.coli
Target Names	SPATA2L
Protein Names	Recommended name: Spermatogenesis-associated protein 2-like protein Short name= SPATA2-like protein
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.