



Recombinant Human WD repeat-containing protein 49 (WDR49)

Product Code	CSB-YP810271HU
Abbreviation	WDR49
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IV35
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAWREKSKKR LNMTSFNIAQ GIHAFDYHSR LNLIATAGIN NKVCLWNPYV VSKPVGVLWG HSASVIAVQF FVERKQLFSF SKDKVLRLLWD IQHQLSIQRI ACSFPKSQDF RCLFHFDEAH GRLFISFNNQ LALLAMKSEA SKRVKSHEKA VTCVLYNSIL KQVISSDTGS TVSFWMIDTG QKIKQFTGCH GNAEISTMAL DANETRLTGT STDGTVKIWD FNGYCHHTLN VGQDGAVDIS QILILKKKIL VTGWERAITV FRPQNFNQFF IQPEEWKGGI QHHDDILCAA FLPPQTLVTG SYDGEIVLWN NSTENAHHVL HPDYQRLLKS KLDTKPQKLL SAGRSQPSHP MADHSTTGVR NFEIDTEGKN AVMLRCFLKA RKNTAVTGGA NLVSCGGSGY VRFWDIYKKQ LLAEFLAHSG VGSIIIMSTDK MNRYLTTGDL DGWLKIWNIE EYCLNSSKNK ITKAPTLIRS FQPHEDRISS LEMCEPGGQL LISSSADCS ICVTGVCNAP VWIFGQAKHW HIENCLFLPK RDTNLVESEI QKEISLFSKE ESCLDPTEHS LLNKKNKDDS TYNVRPSEDI NLDIKYKERS TCMKETQKPY YGEVIKKSFS TFRSLNIGAL EELPEVNKPA FLLDPEKYFR KEPEEERPQI LEAPSLFKTL KAVFDEKNLF PKEILHHERK AKQLCQEKSC EVKKNKK
Source	Yeast
Target Names	WDR49
Protein Names	Recommended name: WD repeat-containing protein 49
Expression Region	1-697
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.