



# Recombinant Human DENN domain-containing protein 1C (DENND1C), partial

<b>Product Code</b>	CSB-EP810272HU-B
<b>Abbreviation</b>	DENND1C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IV53
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESRAEGGSP AVFDWFFEAA CPASLQEDPP ILRQFPPDFR DQEAMQMVPK FCFPFDVERE PPSPAVQHFT FALTDLAGNR RFGFCRLRAG TQSCLCILSH LPWFEVFKL LNTVGDLLAQ DQVTEAEELL QNLFQQSLSG PQASVGLELG SGVTVSSGQG IPPPTRGNSK PLSCFVAPDS GRLPSIPENR NLTELVVAVT DENIVGLFAA LLAERRVLLT ASKLSTLTSC VHASCALLYP MRWEHVLPIPT LPPHLLDYCC APMPYLIGVH ASLAERVREK ALEDVVVLNV DANTLETTFN DVQALPPDVV SLLRLRLRKV ALAPGEGVSR LFLKAQALLF GGYRDALVCS PGQPVTFFSEE VFLAQKPGAP LQAFHRRAVH LQLFKQFIEA RLEKLNKGGEG FSDQFEQEIT GCGASSGALR SYQLWADNLK KGGGALLHSV KAKTQPAVKN MYRSAKSLK GVQSLLMYKD GDSVLQRGGS LRAPALPSRS DRLQQRLPIT QHFGKNRPLR PSRRRQLEEG TSEPPGAGTP PLSPEDGCP WAEALDSSF LGSGEELDLL SEILDSLMSG AKSAGSLRPS QSLDCCHRGD LDSCFSLPNI PRWQPDDKKL PEPEPQPLSL PSLQNASSLD ATSSSKDSRS QLIPSESDQE VTSPSQSSTA SADPSIWGDP KPSPLTEPLI LHLTPSHKAA EDSTAQENPT PWLSTAPTEP SPPEPQILA PTKPNFDIAW TSQPLDPSSD PSSLEDPRAR PPKALLAERA HLQPREPGA LNSPATPTSN CQKSQPSSRP RVADLKKCFE G
<b>Source</b>	E.coli
<b>Target Names</b>	DENND1C
<b>Protein Names</b>	Recommended name: DENN domain-containing protein 1C Alternative name(s): Connecdenn 3 Protein FAM31C
<b>Expression Region</b>	1-801
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.