



Recombinant Human DENN domain-containing protein 1C (DENND1C), partial

Product Code	CSB-MP810272HU
Abbreviation	DENND1C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IV53
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESRAEGGSP AVFDWFFEAA CPASLQEDPP ILRQFPPDFR DQEAMQMVPK FCFPFDVERE PPSPAVQHFT FALTDLAGNR RFGFCRLRAG TQSCLCILSH LPWFEVFKL LNTVGDLLAQ DQVTEAEELL QNLFQQSLSG PQASVGLELG SGVTVSSGQG IPPPTRGNSK PLSCFVAPDS GRLPSIPENR NLTELVVAVT DENIVGLFAA LLAERRVLLT ASKLSTLTSC VHASCALLYP MRWEHVLPIPT LPPHLLDYCC APMPYLIGVH ASLAERVREK ALEDVVVLNV DANTLETTFN DVQALPPDVV SLLRLRLRKV ALAPGEGVSR LFLKAQALLF GGYRDALVCS PGQPVTFFSEE VFLAQKPGAP LQAFHRRAVH LQLFKQFIEA RLEKLNKGEK FSDQFEQEIT GCGASSGALR SYQLWADNLK KGGGALLHSV KAKTQPAVKN MYRSAKSLK GVQSLLMYKD GDSVLQRGGS LRAPALPSRS DRLQQRLPIT QHFGKNRPLR PSRRRQLEEG TSEPPGAGTP PLSPEDGCP WAEALDSSF LGSGEELDLL SEILDSLMSG AKSAGSLRPS QSLDCCHRGD LDSCFSLPNI PRWQPDDKKL PEPEPQPLSL PSLQNASSLD ATSSSKDSRS QLIPSESDQE VTSPSQSSTA SADPSIWGDP KPSPLTEPLI LHLTPSHKAA EDSTAQENPT PWLSTAPTEP SPPEPQILA PTKPNFDIAW TSQPLDPSSD PSSLEDPRAR PPKALLAERA HLQPREPGA LNSPATPTSN CQKSQPSSRP RVADLKKCFE G
Source	Mammalian cell
Target Names	DENND1C
Protein Names	Recommended name: DENN domain-containing protein 1C Alternative name(s): Connecdenn 3 Protein FAM31C
Expression Region	1-801
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.