



# Recombinant Human Protein FAM71A (FAM71A)

<b>Product Code</b>	CSB-BP810307HU
<b>Abbreviation</b>	FAM71A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYT1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNADFLLPYY TAQSGSSMSM FNTTMGKLQR QLYKGEYDIF KYAPIFESDF IQITKRGEVI DVHNRVRMVT MGIARTSPIL PLPDVMLLAR PATGCEEYAG HGQATKRKKR KAAKNLELTR LLPLRFVRIS VQDHEKQQLR LKFATGRSCY LQLCPALDTR DDLFAYWEKL IYLLRPPMES NSSTCGIPAE DMMWMPVFQE DRRSLGAVNL QGKGDQDQVS IQSLHMVSEV CGATSAAYAG GEGQLQDNFNK PTNVLNASIP KTSTELAEPP ATGGIKEAAA AGAAAGAATG TVAGALSVAA ANSAPGQVSA AIAGAATIGA GGNKGNMALA GTASMAPNST KVAVAGAAGK SSEHVSSASM SLSREGSVSL AIAGVVLTSR TAAEADMDDAA AGPPVSTRQS KSSLGQHGGR ERTQASAEGC KEGRERREKD RALGRSSHRR RTGESRHKTR GDKIAQKSSS RSSFSHRANR DDKKEKGCN PGSSRHRDSH KGVSHTPISK ESRTSHKSGR SLWTTSSGSS KGLGRVSSFL RNVRANLTTK VVGTPHGRDV NVMKMAERS TNVAIAETAEGGQGLETVGS MTPDIMETVT FEAH
<b>Source</b>	Baculovirus
<b>Target Names</b>	FAM71A
<b>Protein Names</b>	Recommended name: Protein FAM71A
<b>Expression Region</b>	1-594
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.