



Recombinant Human Protein FAM71A (FAM71A)

Product Code	CSB-EP810307HU-B
Abbreviation	FAM71A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNADFLLPYY TAQSGSSMSM FNTTMGKLQR QLYKGEYDIF KYAPIFESDF IQITKRGEVI DVHNRVRMVT MGIARTSPIL PLPDVMLLAR PATGCEEYAG HGQATKRKKR KAAKNLELTR LLPLRFVRIS VQDHEKQQLR LKFATGRSCY LQLCPALDTR DDLFAYWEKL IYLLRPPMES NSSTCGIPAE DMMWMPVFQE DRRSLGAVNL QGKGDQDQVS IQSLHMVSEV CGATSAAYAG GEGQLQDNFNK PTNVLNASIP KTSTELAEPP ATGGIKEAAA AGAAAGAATG TVAGALSVAA ANSAPGQVSA AIAGAATIGA GGNKGNMALA GTASMAPNST KVAVAGAAGK SSEHVSSASM SLSREGSVSL AIAGVVLTSR TAAEADMDDAA AGPPVSTRQS KSSLGQHGGR ERTQASAEGC KEGRERREKD RALGRSSHRR RTGESRHKTR GDKIAQKSSS RSSFSHRANR DDKKEKGCN PGSSRHRDSH KGVSHTPISK ESRTSHKSGR SLWTTSSGSS KGLGRVSSFL RNVRANLTTK VVGTPHGRDV NVMKMAERS TNVAIAETAEGGQGLETVGS MTPDIMETVT FEAH
Source	E.coli
Target Names	FAM71A
Protein Names	Recommended name: Protein FAM71A
Expression Region	1-594
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.