



Recombinant Human Tensin-4 (TNS4)

Product Code	CSB-MP810319HU
Abbreviation	TNS4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZW8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AP CDEPRRTLHP APSPSLPPQC SYTTTEGWGA QALMAPVPCM GPPGRLQQAP QVEAKATCFL PSPGEKALGT PEDLDSYIDF SLESLNQMIL ELDPTFQLLP PGTGGSQAEI AQSTMSMRKK EESEALDIKY IEVTSARSRC HDGPQHCSSP SVTPPFGSLR SGGLLLSRDV PRETRSSSES LIFSGNQGRG HQRPLPPSEG LSPRPPNSPS ISIPCMGSKA SSPHGLGSPL VASPRLEKRL GGLAPQRGSR ISVLSASPVS DVSYMFGSSQ SLLHSSNSSH QSSSRSLSP ANSSSSLHSL GSVSLCTRPS DFQAPRNPTL TMGQPRTPHS PPLAKEHASS CPPSITNSMV DIPIVLINGC PEPGSSPPQR TPGHQNSVQP GAASPSNPCP ATRSNSQTLS DAPFTTCPEG PARDMQPTMK FVMDTSKYWF KPNITREQAI ELLRKEEPGA FVIRDSSSYR GSFGALKVQ EVPASQSRP GEDSNDLIRH FLIESSAKGV HLGGADEEY FGSLSAFVCQ HSIMALALPC KLTIPQRELG GADGASDSTD SPASCQKKSA GCHTLYLSSV SVETLTGALA VQKAISTTFE RDILPTPTVV HFKVTEQGIT LTDVQRKVFF RRHYPLTTLR FCGMDPEQRK WQKYCKPSWI FGFVAKSQTE PQENVCHLFA EYDMVQPASQ VIGLVTALLQ DAERM
Source	Mammalian cell
Target Names	TNS4
Protein Names	Recommended name: Tensin-4 Alternative name(s): C-terminal tensin-like protein
Expression Region	19-715
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.