



Recombinant Reston ebolavirus Matrix protein VP40 (VP40)

Product Code	CSB-BP810348RCJ
Abbreviation	VP40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8JPX9
Product Type	Recombinant Protein
Immunogen Species	Reston ebolavirus (strain Reston-89) (REBOV) (Reston Ebola virus)
Purity	≥85% (SDS-PAGE)
Sequence	MRRGVLPTAP PAYNDIAYPM SILPTRPSVI VNETKSDVLA VPGADVPSNS MRPVADDNID HSSHTPSGVA SAFILEATVN VISGTKVLMK QIPIWLPLGV ADQKIYSFDS TTAAIMLASV TVTHFGKISN PLVRVNRLGP GIPDHPLRLL RLGNQAFLQE FVLPPVQLPQ YFTFDLTALK LITQPLPAAT WTDETPAGAV NALRPGLSLH PKLRPILLPG KTGKKGHASD LTSPDKIQT MNAIPDLKIV PIDPTKNIVG IEVPELLVQR LTGKKPQPKN GQPIIPVLLP KYVGLDPISP GDLMVITQD CDSCHSPASH PYHMDKQNSY Q
Source	Baculovirus
Target Names	VP40
Protein Names	Recommended name: Matrix protein VP40 Alternative name(s): Membrane-associated protein VP40
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.