



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein SKIP11 (SKIP11)

Product Code	CSB-EP810536DOA-B
Abbreviation	SKIP11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L736
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLED RSPDSC LSTRVFSSSR LSESNWSNSY MYPEDDDKLL GNGKRALEVV GEVRQTKSLK LMGFSIIYDS DSSDYSLSGG EEQADAAIGD GSSSRQEQEQ QSD FNDNGGD SSDSHSLINE IGRDNSIDCL IRCSRSYDYG IASLNRNFRS LVKSGEIYRL RRQNGFVEHW VYFSCQLLEW VAFDPVERRW MQLPTMPSSV TFMCADKESL AVGTDLLVLG KDDFSSHVIY RYSLLTNSWS SGMKMNSPRC LFGSASLGEI AIFAGGCDSQ GKILDFAEMY NSELQTWITL PRMNKPRKMC SGVFM DGK FY VIGGIGGADS KGLTCGEEYD LETKKWTQIP DLSPPRSRAD QADMSPAAEA PPLVAVVNNQ LYAADHADME VRKYDKENKK WLTVGR LPER AGSVNGWGLA FRACGERLIV IGGPKCSGGG FIELNSWIPS DGGPPQW TLL DRKHSPTFVY NCAVMGC
Source	E.coli
Target Names	SKIP11
Protein Names	Recommended name: F-box/kelch-repeat protein SKIP11 Alternative name(s): SKP1-interacting partner 11
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.