



Recombinant Arabidopsis thaliana ATP-dependent DNA helicase Q-like 4A (RECQL4A), partial

Product Code	CSB-EP810545DOA
Abbreviation	RECQL4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L840
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Source	E.coli
Target Names	RECQL4A
Protein Names	Recommended name: ATP-dependent DNA helicase Q-like 4A EC= 3.6.4.12 Alternative name(s): RecQ-like protein 4A Short name= AtRecQ4A Short name= AtRecQl4A SGS1-like protein Short name= AtSGS1
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.