



Recombinant Arabidopsis thaliana Putative clathrin assembly protein At4g40080 (At4g40080)

Product Code	CSB-MP810549DOA
Abbreviation	At4g40080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L936
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MGRITSFADL IGRIKDKASQ SKAALVSSNT KSKTLSFHLS VLRATTHDPS TPPGNRHLAV ILSAGTGSRA TASSAVESIM ERLHTTGDAC VALKSLIIH HIVKHGRFIL QDQLSVFPAS GGRNYLKLSA FRDEKSPLMW ELSSWVRWYA LYLEHLLSTS RIMGFFISST SSTITHKEEYE EMVSSLTNSD LLREIDALVG LLEEACKIPD LPFSGGKSLA DKITQLVGED YVSSINELYT RFNEFKERSN TLSFGDTIEL VCALKRLESC KERLSEICHG NWKRGWIDGF WGLVLEVKGI IGNLEDNYGQ IEKSIVGFGK RDKGYESARF TDRLIIGYSN PVRFSSGRFS NVDRFNFPVS GRVLC
Source	Mammalian cell
Target Names	At4g40080
Protein Names	Recommended name: Putative clathrin assembly protein At4g40080
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.