



Recombinant Thermosynechococcus elongatus Elongation factor Tu (tuf)

Product Code	CSB-MP810877TDQ
Abbreviation	tuf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8DI42
Product Type	Recombinant Protein
Immunogen Species	Thermosynechococcus elongatus (strain BP-1)
Purity	>85% (SDS-PAGE)
Sequence	MARAKFERTK PHVNVGTIGH VDHGKTTLTA AITMTLAAQG KAQARKYDEI DAAPEEKARG ITINTAHVEY ETEKRHYAHV DCPGHADYVK NMITGAAQMD GAILVVAATD GAMPQTKEHI LLARQVGVPS IVVFLNKVDM VDDEELLELV ELELRELLNE YEFPGDEVPI IRGSLKALE AMTANPKTLR GENEWVDKIY ELMDAVDNYI PTPERDVDKP FLMAVEDVFS ITGRGT VATG RIERGRIKLN ETVELVGLRE TRTTT VTGIE MFKKSLEEGI AGDNAGLLLR GLKKEDVERG MVIAKPGSIT PHTKFEGEVY VLTEKEGGRK TPFFAGYRPQ FYVRTTDV TG TITSFTSDDG SAAEMVMPGD RIKMTVELIQ PIAIEQGMRF AIREGGRTIG AGVVSIIIE
Source	Mammalian cell
Target Names	tuf
Protein Names	Recommended name: Elongation factor Tu Short name= EF-Tu
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.