



Recombinant Putative ribosome biogenesis GTPase RsgA (rsgA)

Product Code	CSB-EP811241EGX
Abbreviation	rsgA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FAL3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSKNKLSKGQ QRRVNANHQR RLKTSKEKPD YDDNLFGEPE EGIVISRFGM HADVESADGD VHRCNIRRTI RSLVTGDRVV WRP GKPA AEG VNVKGIVEAV HERTSVLTRP DFYDGVKPIA ANIDQIVIVS AILPELSLNI IDRYLVACET LQIEPIIVLN KIDLLDDEGM EfvNEQMDIY RNIGYRVLMV SSHTQDGLKP LEEALTGRIS IFAGQSGVGK SLLNALLGL QKEILTNDVS DNSGLGQHTT TAARLYHFPH GGDVIDSPGV REFGLWHLEP EQITQGFVEF HDYLG LCKYR DCKHDTDPGC AIREAVEEGK IAETRFENYH RILESMAQVK TRKNFSDTDD
Source	E.coli
Target Names	rsgA
Protein Names	Recommended name: Putative ribosome biogenesis GTPase RsgA EC= 3.6.1.-
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.