



# Recombinant 3-ketoacyl-CoA thiolase (fadA)

<b>Product Code</b>	CSB-EP811249EGX-B
<b>Abbreviation</b>	fadA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FBI3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEQVVIVDAI RTPMGRSKGG AFRNVRAEDL SAHLMRSLLA RNPALAAAL DDIYWGCVQQ TLEQGFNIAR NAALLAEVPH SVPVTVNRL CGSSMQALHD AARMIMTGDA QAQLVGGVEH MGHVPM SHGV DFHPGLSRNV AKAAGMMGLT AEMLARMHGI SREMQDAFAA RSHARAWAAT QSGAFKNEII PTGGHDADGV LKQFNIDEVI RPETTVEALA TLRPAFDPVS GTVTAGTSSA LSDGAAAMLV MSESARDLG LKPRACVRSM AVVGCDPSIM GYGPVPASKL ALKKAGLSAS DIGVFEMNEA FAAQILPCIK DLGLMEQIDE KINLNGGAIA LGHPLGCSGA RISTTLLNLM ERKDVQFGLA TMCIGLGQGI ATVFERV
<b>Source</b>	E.coli
<b>Target Names</b>	fadA
<b>Protein Names</b>	Recommended name: 3-ketoacyl-CoA thiolase EC= 2.3.1.16 Alternative name(s): Acetyl-CoA acyltransferase Beta-ketothiolase Fatty acid oxidation complex subunit beta
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.