



Recombinant Protein CsiD (csiD)

Product Code	CSB-EP811269EGX-B
Abbreviation	csiD
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FES6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MNALTAVQNN AVDSGQDYSG FTLIPSAQSP RLELTFTEQ TTKQFLEQVA EWPVQALEYK SFLRFRVGKI LDDLCANQLQ PLLLKTLLNR AEGALLINAV GIDDVAQADE MVKLATAVAH LIGRSNFDAM SGQYYARFVV KNVDNSDSYL RQPHRVMELH NDGTVEEIT DYVLMMKIDE QNMQGGNSLL LHLDDWEHLD HFFRHPLARR PMRFAAPPSK NVSKDVFHPV FDVDQQGRP V MRYIDQFVQP KDFEEGVWLS ELSDAIETSK GILSVPVPG KLLINLFW LHGRDRFTPH PDLRRELMRQ RGYFAYATHH YQTHQ
Source	E.coli
Target Names	csiD
Protein Names	Recommended name: Protein CsiD
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>