



Recombinant Brucella suis biovar 1 Putative membrane protein insertion efficiency factor (BR1073, BS1330_I1069)

Product Code	CSB-EP811398BMQ-B
Abbreviation	BR1073, BS1330_I1069
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8G0L6
Product Type	Recombinant Protein
Immunogen Species	Brucella suis biovar 1 (strain 1330)
Purity	>85% (SDS-PAGE)
Sequence	MGSCGGKHTG KGAPKPYSRN FTDPWRKTPG RLFGTALIRF YQITLSSLIG NSCRHLPTCS EYAYEAIARH GLWRGGWMGF FRVVRGPFPG THGFDPVPRE LSPDLKWYMP WRYWRCSASR TGK
Source	E.coli
Target Names	BR1073
Protein Names	Recommended name: Putative membrane protein insertion efficiency factor
Expression Region	1-123
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.