



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 29

(Os03g0207400, LOC_Os03g10950)

| | |
|--------------------------|--|
| Product Code | CSB-EP811496OFG-B |
| Abbreviation | Os03g0207400, LOC_Os03g10950 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8H063 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MGALRRWLPC CCCCCRGGGG GGGGGSVGDG LVWDVALKAH ASGDYSVAVA QANEALDQA QVVFSPAATL VGVYDGHGGP EAARFVNKRL FSLIQEFAAQ SGGISAEVLE KAFGETEEEF VASVQRSWPS QPRILSVGSC CLVGAIEDGT LYVANLGDGR AVLGRRSAAG AAHGRKGKNR VVPERLSRDH NVADEDVRRE LKELHPDDSH IVLNTHGVWR IKGIIQVSRS IGDVYLKKPE ICKSNPMLQQ TICPFPLRRP VMSAVPTIKT RKL RPDGQFV IFASDGLWEQ LTDEAAVAIV AGSPRRGVAM RLVRAAQLEA ARKKDKKYER IRTIEKGQRR HFHDDITVVV LFLDKCRGKA GRGDEIDGTD GPVDVFLSP DDREDPTRPV LR |
| Source | E.coli |
| Target Names | Os03g0207400 |
| Protein Names | Recommended name: Probable protein phosphatase 2C 29 Short name= OsPP2C29 EC= 3.1.3.16 |
| Expression Region | 1-392 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.