



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 29 (Os03g0207400, LOC_Os03g10950)

Product Code	CSB-MP811496OFG
Abbreviation	Os03g0207400, LOC_Os03g10950
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H063
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGALRRWLPC CCCCCRGGGG GGGGGSVGDG LVWDVALKAH ASGDYSVAVA QANEALDQA QVVFSPAATL VGVYDGHGGP EAARFVNKRL FSLIQEFAAQ SGGISAEVLE KAFGETEEEF VASVQRSWPS QPRILSVGSC CLVGAIEDGT LYVANLGD SR AVLGRRSAAG AAHGRKGKNR VVPERLSRDH NVADEDVRRE LKELHPDDSH IVLNTHGVWR IKGIIQVSRS IGDVYLKKPE ICKSNPMLQQ TICPFPLRRP VMSAVPTIKT RKL R PGDQFV IFASDGLWEQ LTDEAAVAIV AGSPRRGVAM RLVRAAQLEA ARKKDKYER IRTIEKGQRR HFHDDITVVV LFLDKCRGKA GRGDEIDGTD GPVDVFLSP DDREDPTRPV LR
Source	Mammalian cell
Target Names	Os03g0207400
Protein Names	Recommended name: Probable protein phosphatase 2C 29 Short name= OsPP2C29 EC= 3.1.3.16
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.