



Recombinant Human TIR domain-containing adapter molecule 1 (TICAM1), partial

Product Code	CSB-EP811598HU1
Abbreviation	TICAM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUC6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RLDEHSQIFARKVANTFKPHRLQARKAMWRKEQDTRALREQSQHLDGERMQ AAALNAAYS AYLQSYLSYQAQMEQLQVAFGSHMSFGTGAPYGARMPFGGQVPLGAPPPFP TWPGCPQPP PLHAWQAGTPPPSPQPAAFPQSLPFPQSPAFPTASPAPPQSPGLQPLIIHHA QMVQLGL NNHMWNQRGSQAPEDKTQEAE
Source	E.coli
Target Names	TICAM1
Protein Names	Recommended name: TIR domain-containing adapter molecule 1 Short name= TICAM-1 Alternative name(s): Proline-rich, vinculin and TIR domain-containing protein B Putative NF-kappa-B-activating protein 502H Toll-interleukin-1 receptor d
Expression Region	512-712
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.