



Recombinant Human SURP and G-patch domain-containing protein 1 (SUGP1)

Product Code	CSB-EP811617HU-B
Abbreviation	SUGP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IWZ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLKMDNRDV AGKANRWFGV APPKSGKMMN NILHQEELIA QKKREIEAKM EQKAKQNQVA SPQPPHPGEI TNAHNSSCIS NKFANDGSFL QQFLKLQKAQ TSTDAPTSAP SAPPSTPTPS AGKRSLISR RTGLGLASLP GPVKSYSYSHAK QLPVAHRPSV FQSPDEDEEE DYEQWLEIKV SPPEGAETRK VIEKLARFVA EGGPELEKVA MEDYKDNPAF AFLHDKNSRE FLYYRKKVAE IRKEAQKSQA ASQKVSPPED EEVKNLAEKL ARFIADGGPE VETIALQNNR ENQAFSFLYE PNSQGYKYYR QKLEEFRKAK ASSTGSFTAP DPGLKRKSP EALSGSLPPA TTCPASSTPA PTIIPAPAAP GKPASAATVK RKRKSRWGPE EDKVELPPAE LVQRDVDASP SPLSVQDLKG LGYEKGKPVG LVGVTELSDA QKKQLKEQQE MQQMYDMIMQ HKRAMQDMQL LWEKAVQQHQ HGYDSDEEVD SELGTWEHQL RRMEMDKTRE WAEQLTKMGR GKHFIDFLP PDELEKFMET FKALKEGREP DYSEYKEFKL TVENIGYQML MKMGWKEGEG LGSEGQGIKN PVNKGTTTVD GAGFGIDRPA ELSKEDDEYE AFRKRMMLAY RFRPNPLNNP RRPYY
Source	E.coli
Target Names	SUGP1
Protein Names	Recommended name: SURP and G-patch domain-containing protein 1 Alternative name(s): RNA-binding protein RBP Splicing factor 4
Expression Region	1-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.