



Recombinant Human Transcription factor Sp8 (SP8)

Product Code	CSB-BP811623HU
Abbreviation	SP8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXZ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLAATCNKIG SPSPSPSSLS DSSSSFGKGF HPWKRSSSSS SASCNVVGSS LSSFGVSGAS RGGSSSSAAA AAAAAAAAAA ALVSDSFSCG GSPGSSAFSL TSSSAAAAAA AAAAAASSSP FANDYSVFQA PGVSGGSGGG GGGGGGGSSA HSQDGSHQPV FISKVHTSVD GLQGIYPRVG MAHPYESWFK PSHPGLGAAG EVGSAGASSW WDVGAGWIDV QNPNSAAALP GSLHPAAGGL QTSLHSPGG YNSDYSGLSH SAFSSGASSH LLSPAGQHLM DGFKPVLPGS YPDSAPSPLA GAGGSMLSAG PSAPLGGSPR SSARRYSGRA TCDCPNCQEA ERLGPAGASL RRGKGLHSCHI PGCGKVIYGKT SHLKAHLRWH TGERPFVCNW LFCGKRFRTR DELQRHLRTH TGEKRFACPV CNKRFMRSDH LSKHVKTHSG GGGGGGSAGS GSGGKKGSDT DSEHSAAGSP PCHSPELLQP PEPGHRNGLE
Source	Baculovirus
Target Names	SP8
Protein Names	Recommended name: Transcription factor Sp8 Alternative name(s): Specificity protein 8
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.