



# Recombinant Human Kelch repeat and BTB domain-containing protein 2 (KBTBD2)

<b>Product Code</b>	CSB-EP811624HU
<b>Abbreviation</b>	KBTBD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IY47
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTQDERQIN TEYAVSLEEQ LKLFYEQQLF TDIVLIVEGT EFPCHKMVLV TCSSYFRAMF MSGLSSESKQT HVHLRNVDAAL TLQIIITYAY TGNLAMNDST VEQLYETACF LQVEDVLQRC REYLIKKINA ENCVRLLSFA DLFSCEELKQ SAKRMVEHKF TAVYHQDAFM QLSHDLLIDI LSSDNLNVEK EETVREAAML WLEYNTESRS QYLSSVLSQI RIDALSEVTQ RAWFQGLPPN DKSVVVQGLY KSMKFFKPR LGMTKEEMMI FIEASSENPC SLYSSVCYSP QAEKVYKLCV PPADLHKVGT VVTPDNDIYI AGGQVPLKNT KTNHKSLSKL QTAFRTVNCV YWFDAQQNTW FPKTPMLFVR IKPSLVCCEG YIYAIGGDSV GGELNRRTVE RYDTEKDEWT MVSPPLCAWQ WSAAVVVHDC IYVMTLNLMY CYFPRSDSWV EMAMRQTSRS FASAAAFGDK IFYIGGLHIA TNSGIRLPSG TVDGSSVTVE IYDVNKNEWK MAANIPAKRY SDPCVRAVVI SNSLCVFMRE THLNERAKYV TYQYDLELDR WSLRQHISER VLWDLGRDFR CTVGKLYPSC LEESPWKPPT YLFSTDGTEE FELDGEMVAL PPV
<b>Source</b>	E.coli
<b>Target Names</b>	KBTBD2
<b>Protein Names</b>	Recommended name: Kelch repeat and BTB domain-containing protein 2 Alternative name(s): BTB and kelch domain-containing protein 1
<b>Expression Region</b>	1-623
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.