



# Recombinant Human Stromal membrane-associated protein 1 (SMAP1)

<b>Product Code</b>	CSB-BP811625HU
<b>Abbreviation</b>	SMAP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYB5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MATRSCREKA QKLNEQHQLI LSKLLREEDN KYCADCEAKG PRWASWNIGV FICIRCAGIH RNLGVHISRV KSVNLDQWTA EQIQCMQDMG NTKARLLYEA NLPENFRRPQ TDQAVEFFIR DKYEKKKYD KNAIAITNIS SSDAPLQPLV SSPSLQAAVD KNKLEKEKEK KKEEKREKE PEKPAKPLTA EKLQKKDQQL EPKKSTSPKK AAAPTVDLLG LDGPAVAVPT NGNTTVPLN DDLDFGPMI SNPLPATVMP PAQGTPSAPA AATLSTVTSG DLDFLTEQTT KSEEVAKKQL SKDSILSLYG TGTIQQQSTP GVFMGPTNIP FTSQAPAAFQ GFPSMGVVPV AAPGLIGNVM GQSPSMMVGM PMPNGFMGNA QTGVMLPQN VVGPPQGGMVG QMGAPQSKFG LPQAQQPQWS LSQMNQQMAG MSISSATPTA GFGQPSSTTA GWSGSSSGQT LSTQLWK
<b>Source</b>	Baculovirus
<b>Target Names</b>	SMAP1
<b>Protein Names</b>	Recommended name: Stromal membrane-associated protein 1
<b>Expression Region</b>	1-467
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is similar to the mouse stromal membrane-associated protein-1. This similarity suggests that this human gene product is also a type II membrane glycoprotein involved in the erythropoietic stimulatory activity of stromal cells. Alternate splicing results in multiple transcript variants encoding different isoforms.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.