



Recombinant Human Ubiquilin-like protein (UBQLNL)

Product Code	CSB-YP811630HU
Abbreviation	UBQLNL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYU4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWHAI S R T S R MSQSGCPSGL LADKNISSA TRVIVKTAGN QKDFMVADDI SVRQFKEMLL AHFQCQMDQL VLVFMGCLLK DHDTLSQ R GI MDGHTIYLVI KSKQGS R SLA H S FRDLPTND PCH R DRNTKG NSSRVHQPTG MNQAPVELAH FVGSDAPKVH TQNLEVSHPE CKAQMLENPS IQRLLSNMEF MWQFISEHLD TQQLMQQNPE VSRLLLDNSE ILLQTLELAR NLAMIQEIMQ IQQPSQNLEY PLNPQPYLGL ETMPGGNNAL GQNYADINDQ MLNSMQDPFG GNPFTALLAG QVLEQVQSSP PPPPPSQEQQ DQLTQHPATR VIYNSSGGFS SNTSANDTLN KVNHTSKANT AMISTKGQSH ICATRQPAWI PALPSIELTQ QLQEEYKDAT VSLSSSRQTL KGD L QLSDEQ SSSQITGGMM QLLMNNPYLA AQIMLFTSMP QLSEQWRQQL PTFLQQTQIS DLLSA
Source	Yeast
Target Names	UBQLNL
Protein Names	Recommended name: Ubiquilin-like protein
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.