



Recombinant Human TBC1 domain family member 21 (TBC1D21)

Product Code	CSB-BP811631HU
Abbreviation	TBC1D21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYX1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTTLSPENSL SARQSASFIL VKRKPPIDKT EWDSFFDESG HLAKSRDFIC VNILERGLHP FVRTEAWKFL TGYFSWQSSQ DERLTVDSMR RKNYKALCQM YEKIQPLEN LHRNFTETRN NIARDIQKIY DKDPLGNVLI DKKRLEKILL LSYVCNTQAE YQQGFHEMMM LFQLMVEHDH ETFWLFQFFL QKTEHSCVIN IGVAKNLDML STLITFLDPV FAEHLKKGKA GAVQSLFPWF CFCFQRAFKS FDDVWRLWEV LLTGKPCRNF QVLVAYSMLQ MVREQVLQES MGGDDILLAC NNLIDLDADE LISAACVVYA ELIQKDVPQT LKDDFFL
Source	Baculovirus
Target Names	TBC1D21
Protein Names	Recommended name: TBC1 domain family member 21
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.