



Recombinant Human Proline-rich protein 22 (PRR22)

Product Code	CSB-MP811633HU
Abbreviation	PRR22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZ63
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQHPKPFCAP AAPQEGFSPQ SLEGAEVLGN QPAPTCAEPP PAMGSLNLYH PPDPEKEVFP APPAGFQMAP CGCFFDPRIY RIEWTTPLDLG QSALYKLAAS SGGPAGVPSA PGSYLLEPQP YLKAPGLPPY PHYQQAPGGP QFLLPYFPPE GPGPEALGFV GDAGPAAFVE LPLPPLEEGP APLPPPPPKE NKPPPVLITL PAEPTLPPDA YSHLQGH LGH FPGPEPLAFP VKELQGSGAR PGVPLYPPGL SELKVAEVKE GALLGAGKAK APKTARALAL PDKVLLEDAM KLFDCPLPGAS EPEGTLCEVP GPALPDSSGG NSADDIRSLC LPEELLSFDY SVPEILDTVS NVDYFFNFKA LDEEQPPHPG PPATNTPAPI LSGKRKASTA KKGKPGRKAR QPAGPASATP GGPREDLGAT PH</p>
Source	Mammalian cell
Target Names	PRR22
Protein Names	Recommended name: Proline-rich protein 22
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.