



Recombinant Rat Protein FAM5C (Fam5c)

Product Code	CSB-MP811686RA
Abbreviation	Fam5c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K1M7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	AVSDQHA TSPFDWLLSD KGPFHRSQEY TDFVDRSRQG FSTRYKIYRE FGRWKVNNLA VERRNFLGSP LPLAPEFFRN IRLLGRRPTL QQITENLIKK YGTHFLLSAT LGGEESLTIF VDKRKLSKRP EGSETSTNSS SVTLETLHQL AASYFIDRDS TLRLHHIQI ASTAIKVTTET RTGPLGCSNY DNLDVSSVL VQSPENKIQL QGLQVLLPDY LQERFVQAAL SYIACNSEGE FICKENDCWC LCGPKFPECN CPSMDIQAME ENLLRITETW KAYNSDFEDS DEFKFFMKRL PMNYFLNTST IMHLWTMDSN FQRRYELEN SMKQLFLKAH RIVHKLFSL KRCHKQPLIS LPRQRTSTYW LTRIQSFLYC NENGLLGSFS EETHSCTCPN DQVVCTAFLP CTVGDASACL TCAPDNRTRC GTCNTGYMLS QGLCKPEVAE STDHYIGFET DLQDLEMKYL LQKTDRRIEV HAIFISNDRM LNSWFDPSWR KRMLLTLKSN KYKSSLVHMI LGLSLQICLT KNSTLEPVLA VYINPFGGSH SESWFMPVSE SSFPDWERTK LDLPLQCYNW TLTLGNKWKT FFETVHIYLR SRIKANGPNS NESIYYEPL FIDPSRNLGY MKINNIQVFG YSMHFDPEAI RDLILQLDYP YTQGSQDSAL LQLLEIRDRV NKLSPPGQRR LDLFSCLLRH RLKLSTSEVV RIQSALQAFN AKLPNTVDYD TTKLCS
Source	Mammalian cell
Target Names	Brinp3
Protein Names	Recommended name: Protein FAM5C Alternative name(s): BMP/retinoic acid-inducible neural-specific protein 3
Expression Region	34-766
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.