



Recombinant Mouse Protein phosphatase Slingshot homolog 3 (Ssh3)

Product Code	CSB-EP811700MO
Abbreviation	Ssh3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K330
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ALVTVSRSP PASGHSTPVG PTQDRVRRR GRLQRRQSFA VLRGAVLGLQ DGGDSNVASE ADSEPMEEPS GEEQPTEDQT DKGQGLQSPW KQVQKRHLHL MVELLRPQDD IRLAAQLEAA RPPRLRYLLV VSTGEELSEE AILLGVDFPD SSSHCTLGL VLPLWSDTQV YLDGDGGFSV TSGGQSRIFK PVSQTMWAT LQVLHQACEV ALGSGLVP GG SALAWATYYQ EKLNSDQGCL NEWMAMSDLE SFRPPNAEPG QASEQEKMEQ AILAEWQVL DTSDLDSVTS KEIRQALELR LGCPLQQYRD FIDNQMLLLM AQQDRASRIF PHLYL GSEWN AANLEELQKN RVSHILNMAR EIDNFFPERF TYYNVRVWDE ESAQLLPHWK ETHRFIEDAR AQGTRVLVHC KMGVSRSAAT VLAYAMKQYG WDLEQUALIHV QELRPIVRPN HGFLRQLRTY QGILTASRQS HVWEQKVG VV SPEEPLAPEV STPLPPLPPE PGGSGEVMVM GLEGSQETPK EELGLRPRIN LRGVMRSISL LEPSESESTP EAGGLPEVFS SDEEPLHPFS QLSRAKGGQR VRKGPWPALK SRQSVVALHS AALVASRTRA FQEQQGQEQ SEPGMSSTPR LRKVMRQASV DDSREEDKA
Source	E.coli
Target Names	Ssh3
Protein Names	Recommended name: Protein phosphatase Slingshot homolog 3 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): SSH-like protein 3 Short name= SSH-3L Short name= mSSH-3L
Expression Region	2-649
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.