



# Recombinant Mouse Protein phosphatase Slingshot homolog 3 (Ssh3)

<b>Product Code</b>	CSB-MP811700MO
<b>Abbreviation</b>	Ssh3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K330
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ALVTVSRSP PASGHSTPVG PTQDRVVRRR GRLQRRQSFA VLRGAVLGLQ DGGDSNVASE ADSEPMEEPS GEEQPTEDQT DKGQGLQSPW KQVQKRHLHL MVELLRPQDD IRLAAQLEAA RPPRLRYLLV VSTGEELSEE AILLGVDFPD SSSHCTLGL VLPLWSDTQV YLDGDGGFSV TSGGQSRIFK PVSQTMWAT LQVLHQACEV ALGSGLVP GG SALAWATYYQ EKLNSDQGCL NEWMAMSDLE SFRPPNAEPG QASEQEKMEQ AILAEWQVL DTSDLDSVTS KEIRQALELR LGCPLQQYRD FIDNQMLLLM AQQDRASRIF PHLYL GSEWN AANLEELQKN RVSHILNMAR EIDNFFPERF TYYNVRVWDE ESAQLLPHWK ETHRFIEDAR AQGTRVLVHC KMGVSRSAAT VLAYAMKQYG WDLEQALIHV QELRPIVRPN HGFLRQLRTY QGILTASRQS HVWEQKVG VV SPEEPLAPEV STPLPPLPPE PGGSGEVMVM GLEGSQETPK EELGLRPRIN LRGVMRSISL LEPSESESTP EAGGLPEVFS SDEEPLHPFS QLSRAKGGQR VRKGPWPALK SRQSVVALHS AALVASRTRA FQEQQGQGEQ SEPGMSSTPR LRKVMRQASV DDSREEDKA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ssh3
<b>Protein Names</b>	Recommended name: Protein phosphatase Slingshot homolog 3 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): SSH-like protein 3 Short name= SSH-3L Short name= mSSH-3L
<b>Expression Region</b>	2-649
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.