



Recombinant Arabidopsis thaliana Branched-chain-amino-acid aminotransferase-like protein 3, chloroplastic (At5g57850)

Product Code	CSB-BP811851DOA
Abbreviation	At5g57850
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L493
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	SSY EVGERLKLAR GGQQFLAMYS SVVDGITTDPAAMVLPLDDH MVHRGHGVFD TALIINGYLY ELDQHLDRI LRSASMAKIPL PFDRETIKRI LIQTVSVSGC RDGSLRYWLS AGPGDFLLSP SQCLKPTLYA IVIKTNFAIN PIGVKVVTSS IPIKPPEFAT VKSVNYLPNV LSQMEAEAKG AYAGIWWCKD GFIAEGPNMN VAFVVNGGKE LVMPRFDNVL SGCTAKRTL LAEQLVSKGI LKT VKVMDVT VEDGKKADEM MLIGSGIPIR PVIQWDEEFI GEGKEGPIAK ALLDLLLED M RSGPPSVRVL VPY
Source	Baculovirus
Target Names	DAAT
Protein Names	Recommended name: Branched-chain-amino-acid aminotransferase-like protein 3, chloroplastic
Expression Region	58-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.