



# Recombinant N-acetylneuraminate epimerase 2 (nanM2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP811982EGX                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Abbreviation</b>      | nanM2                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                       |
| <b>Uniprot No.</b>       | Q8CVG5                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Sequence</b>          | S VLPETPVPEFK SGTGVIDNDT VYIGLGSAGT AWYKLDQAK DKRWTALAAF<br>PGGPRDQATS AFIDGNLYVF GGIGKNSEGL TQVFNDVHKY NPKTNSWVKL<br>MSHAPMG MAG HVT FVHNGKA YVTGGVNQNI FNGYFEDLNE AGKDSTAIK<br>INAYYFDKKA EDYFFNKFL SFD PSTQQWS YAGESPWYGT AGAAVVNKG D<br>KTW LINGEAK PGLRTDAVFE LDFTGN NLKW NKLAPVASPD GVAGGFAGMS<br>NDSLIFAGGA GFKGSRENYQ NGKNYAHEGL KKSYSADIHL WHNGKWDKSG<br>ELS QGRAYGV SLLWNN SLLI IGGEAAGGKA VTDSVLISVK DNKVTVQN |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | nanM2                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Names</b>     | Recommended name: N-acetylneuraminate epimerase 2 EC= 5.1.3.24<br>Alternative name(s): N-acetylneuraminate mutarotase 2 Short name= Neu5Ac mutarotase 2 Sialic acid epimerase 2                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 20-368                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                      |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.              |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                       |