



Recombinant Human Heterogeneous nuclear ribonucleoprotein H (HNRNPH1)

Product Code	CSB-EP010609HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31943
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLGTEGGEG FVVKVRGLPW SCSADEVQRF FSDCKIQNGA QGIRFIYTRE GRPSGEAFVE LESEDEVKLA LKKDRETMGH RYVEVFKSNN VEMDWVLKHT GPNSPDTAND GFVRLRGLPF GCSKEEIVQF FSGLEIVPNG ITLPVDFQGR STGEAFVQFA SQEIAEKALK KHKERIGHRY IEIFKSSRAE VRTHYDPPRK LMAMQRPGPY DRPGAGRGYN SIGRGAGFER MRRGAYGGGY GGYDDYNGYN DGYGFGSDF GRDLNYCFSG MSDHRYGDGG STFQSTTGHC VHMRLPYRA TENDIYNFFS PLNPVRVHIE IGPDGRVTGE ADVEFATHED AVAAMSKDKA NMQHRYVELF LNSTAGASGG AYEHRVELF LNSTAGASGG AYGSMQMMGGM GLSNQSSYGG PASQQLSGGY GGYGGQSSM SGYDQVLQEN SSDLFQSNIA
Source	E.coli
Target Names	HNRNPH1
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein H Short name= hnRNP H Cleaved into the following chain: 1. Heterogeneous nuclear ribonucleoprotein H, N-terminally processed
Expression Region	2-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.