



Recombinant *Leptospira interrogans* serogroup Icterohaemorrhagiae serovar Lai 50S ribosomal protein L35 (rpml)

Product Code	CSB-BP812504LOS
Abbreviation	rpml
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8F6Q8
Product Type	Recombinant Protein
Immunogen Species	<i>Leptospira interrogans</i> serogroup Icterohaemorrhagiae serovar Lai (strain 56601)
Purity	>85% (SDS-PAGE)
Sequence	MPKLLKTNRAA AKRFKFTKNN KIKRKS MNTR HILTKKGP KR RRRLRGLTLV HNSDWKSIVR LMPYGV R
Source	Baculovirus
Target Names	rpml
Protein Names	Recommended name: 50S ribosomal protein L35
Expression Region	1-67
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.